## Help sheet 21 - stacks and piles



## Using stacks or piles.

This is done as above by creating a Cartesian ROW. The dimensions of the stack might not be so well known.

Normally you would have some sort of holder to restrain the items in as on the right. However if you are placing objects on a pile with no requirement to pick them up again then a loose pile will work.

Suppose you have a stack that must hold up to 20 plastic microplates like those on the right.

- 1 First stack 20 plates in the holder.
- 2 Create a route 20 lines high reserve at least 21 in case you want approach position.
- 3 Using Jog take the robot to the top plate and learn line 20
- 4 Move robot away and remove all but the bottom plate.
- 5 Using Jog teach line 1.
- 6 Click interpolate.
- 7 If you want an approach position, for example to carefully lower a plate on to the one below it then create the approach as above, e.g.

0 0 100 MOVE for 10mm

then click set approach.

X and Y,

an

```
You can then unload the stack very simply like this.
: FILLSTACK
                 ( STACK is the name of the route
STACK
ALIGN
20 1 DO
                 ( whatever that may entail
    GETPLATE
                 ( go to the tray, next position
    I INTO
                 ( put plate in the tray
    UNGRIP
    UP
L<sub>00</sub>P
Or empty a stack like this:
: GETSTACK
                 ( STACK is the name of the route
STACK
ALIGN
20 1 DO
    I INTO
                 ( go to the tray
    GRIP
    UP
                 ( whatever that may entail
    PUTPLATE
L<sub>0</sub>OP
```

## Help sheet 21 – stacks and piles



Of course you may have more than one stack. In that case you would need a variable for each stack, to keep count.

```
USER STK1
USER STK2
: GET1
           ( name of route for stack 1 - don't forget only the first 5
STAK1
characters are recognized. STACK1 and STACK2 would not work. STK1 @ ( get the count for stack1 from STK1
           ( go to the position noted in STK1
INTO
GRIP
UP
-1 STK1 +! ( decrement the value in STK1
then PUTPLATE or whatever is next.
: GET2
           ( name of route for stack 1 - don't forget only the first 5
STAK2
characters are recognized. STACK1 and STACK2 would not work. STK2 @ ( get the count for stack1 from STK1
           ( go to the position noted in STK1
INTO
GRIP
UP
-1 STK2 +! ( decrement the value in STK2
CR ." fill stacks with plates. Press enter. " KEY DROP
20 STK1 ! 20 STK2 !
Or to put a plate back:
: PUT1
           ( name of route for stack 1 - don't forget only the first 5
STAK1
characters are recognized. STACK1 and STACK2 would not work.
           ( get the count for stack1 from STK1
STK1 @
INTO
           ( go to the position noted in STK1
UNGRIP
UP
1 STK1 +! ( decrement the value in STK1
Make sure you don't over-fill the stack
: PUT1
           ( name of route for stack 1 - don't forget only the first 5
STAK1
characters are recognized. STACK1 and STACK2 would not work.
          ( get the count for stack1 from STK1
DUP 19 > IF 101 ABORT THEN
           ( go to the position noted in STK1
INTO
UNGRIP
UP
1 STK1 +! ( decrement the value in STK1
If your program or user tries to put in more than 20 the system will abort. Enter ERR ? and you will see
101 (or whatever code you allocate) which tells you why the system stopped.
```